ARTZT
neuro

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What is neuro training?

Body and mind are a single entity. Our brain regulates and controls everything in the body. It processes external and internal information (input) with the existing knowledge and creates an end result (output). Particularly in the field of sport, the key output is movement, which depends on the processes taking place in the brain and nervous system. The goal of neuro athletics is to improve performance and prevent injury by optimising these processes. The exercises presented here also enable general health improvements to be achieved. Since all exercises focus on the processes that take place in the brain and nervous system, this form of training is also referred to as neurocentric training.

When developing our tools, we collaborate closely with the internationally renowned expert in neurocentric rehabilitation and athletic performance, Eric Cobb (Z-Health Performance Solutions).

Neuro-training tools

Our neuro-training set consists of visual charts (VOR, peripheral vision, near-far vision, saccades), Brock String, 2 Vision Sticks, Marsden Ball, Pinhole Glasses and an Eye Patch.

Simple and versatile

Since they are so easy to handle, the neuro-training tools are suitable for use in sports, everyday life, work and therapy.

Transportable

Neuro-training tools are small and lightweight and can fit into any bag. This makes them quick and easy to use in any imaginable setting.

Effective

In the case of frequently unrecognised problems in the field of visual perception, balance or proprioception, effective results can be achieved with small daily exercises.

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SET INCLUDES | CONTENU DU SET | CONTENUTO DEL SET

ARTZT
neuro

Vision Charts – train the eye muscles and peripheral perception. (VOR Chart, „Peripheral Awareness“ Chart, Near Far Charts, Saccade Charts)

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Eye Patch – can increase the level of difficulty of the exercises and improve vision limitations by the exercises being done with the strong eye covered up.

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Pinhole Glasses – Relax the eye muscles when focussing. (Please read restrictions of use on packaging!)

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Brock String – Trains synchronous, binocular vision.

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Vision Sticks – Train the targeted control of fixation points as well as the contraction and relaxation of the eye muscles.

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Marsden Ball – trains the eye tracking movements and thus constant focussing on a moving object.

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Neuro-training tools can be used for sports, e-sports, therapy, rehabilitation, fitness, everyday life and at the office. The tools have a multitude of possible uses:

- Optimisation of ocular functioning / visual perception
- Optimisation of vestibular processes and, among other things, orientation in space
- General improvement of the sensory system
- Improvement of orientation and muscular coordination in children
- Training of eye-hand coordination
- Improvement of responsiveness
- Relaxation of overworked structures
- Reduction of burnout, overwork and stress symptoms
- Prevention of injury through optimised information processing
- Improved control after injuries

Complex mental and physical demands

In sport, not only do we move our bodies, we constantly have to check our position in space, as well as that of any opponents nearby, and coordinate them with our own position and, ideally, respond to unforeseeable events, such as uneven ground or collisions with opponents, as best as possible.

The smoother and more precise these processes are, the better the result.

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Example: Fast break in handball

1. Requirements for the visual system (information about the environment):
 - a. Constant switching between close-up and distance vision with the fastest possible accommodation of the eyes, so that the ball, teammates, opponents and the goal can be seen clearly in rapid succession.
 - b. Optimal binocular vision and stable jumps in gaze to fix the ball, goal and opponents and to perceive them exactly where they are at the moment.
 - c. Peripheral field of vision: the bigger this is, the better my peripheral perception is, the faster I can spot upcoming dangers (such as tackles by an opponent).
2. Requirements for the vestibular system (information about position and movement in space):
 - a. Where am I in the room? Where is up and down, for example, during and after a jump shot?
 - b. Rapid direction detection and movement adaptation with appropriate fixation of the eye movement on the respective visual target.
 - c. Realignment, stabilisation or straightening of the body during tackles by an opponent.
 - d. Estimation of the distance to the opponent, teammate or goal.
3. Requirements for the proprioceptive system (self perception):
 - a. The best possible, i.e. precise perception of joint position, muscle tension, speed of movement, tension and pressure are used for rapid error correction and for coordinated movement that is optimally aligned to the situation.
 - b. Maintaining balance even when the opponent tackles.
 - c. Correct movement even without or with limited visual support.

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The exercises presented here serve to prevent injuries and improve performance in sports by processing information as fast as possible. Moreover, the exercises can also lead to improvements in general health beyond sport.

The manual focuses on exercises with small devices for the visual and balance systems. It makes sense to complement comprehensive neurocentric training with exercises to improve proprioception.

Visual system

Almost all of our movements are controlled through the eyes or the visual system. All the important information about our environment is communicated to our brain via the input from the eyes in order to achieve purposeful movement. Even minor information losses or faults can lead to a less than optimal motion result. In addition, strained eyes can lead to headaches, difficulty concentrating and neck tension. They are also closely related to the sense of balance.

Visual training can be thought of as a kind of fitness and coordination training for the eyes. Regular training of this aspect can lead to impressive results. As with any training, however, breaks should also be taken. Visual training is strenuous at first. 10 minutes a day is therefore enough at the beginning. The training can be made easier or harder with variations:

1. Easier: Larger letters or numbers / smaller range of motion / shorter duration / slower work. Rhythm or repeating aloud what has been read can help
2. More difficult: Smaller letters or numbers / greater range of motion / longer duration / speed increase / one-eyed (eyepatch) work / add movement / work with moving objects.

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1. Basic training

Hold the VOR Chart (VOR = vestibulo-ocular reflex) at eye level or hang it on the wall immediately in front of you. With your head remaining absolutely still, trace the lines on the board to the end with your gaze. Hold to the right, hold for 5-10 seconds, to the left, hold, up, hold, down, hold, diagonally to the top right, and so on.

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
2. Peripheral perception / lateral field of vision

Neutral, relaxed state. Hold the Chart „Peripheral Awareness“ about 30-40 cm in front of your face, centred on it. Fix the small „z“ in the middle of the chart. Now read the letters from the inside to the outside, without moving the view away from the centre. Try the exercise one-eyed with the Eye Patch. 2. Peripheral perception / lateral field of vision

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Nah-Tafel

Near Chart
Tableau de proximité
Pannello vicino

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3. Accommodation / focussing ability

Hang the Far Chart against the wall and move so far that you can still see the letters clearly. Likewise, hold the Near Chart at a distance where you can still see the letters clearly. Now read near and far letters alternately. Do not change chart until you have fully recognised the letter. Variations: With Eye Patch / one-eyed

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4. Focussing / stable eye jumps (Saccades)

Hang the two Saccade Charts next to each other at eye level on the wall at a height of approx. 30-40 cm and stand in front of it in the middle (wall distance approx. 1 m). Now jump from letter to letter with your eyes (1st line left, 1st line right, 2nd line left, 2nd line right, etc.). Do not jump on until you can clearly see the letter. Reading the letters aloud supports the exercise. Variations: This exercise can also be performed using Vision Sticks.

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
5. Eye tracking movements / Controlled object tracking

Lie on you back on the floor and relax. Hang the Marsden Ball centrally about 60-100 cm above your face. Swing the ball in all directions and follow it with your eyes (fixating on the letter on the bottom of the ball). Variations: The exercise can also be performed in different stance positions and with different ball heights.

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1. Combination basic exercises (eye tracking movement + final fixation)

Affix the Eye Patch to one eye and throw the Marsden Ball and catch it. The task is to actively track the ball with the open eye, i.e. also until the final catch in the hand with a view of the ball / letters in the palm of the hand. Sticks.

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Vision Stick

Vision Stick
Vision Stick
Vision Stick

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2. Flexing and relaxing

Hold the Vision Stick in front of your face at a distance of about 30-40 cm. Keep your head straight and fixate on a letter on the stick. Guide the stick slowly towards your nose and try to keep the image sharp. Once you see the stick in double or it is very blurry, reverse the movement. If the two eyes are not in sync, the deviant eye can be trained separately with the aid of the eye patch.

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
3. Two-eyed, synchronous vision (binocular)

Hold one end of the Brock String at the tip of the nose and tighten the string. Start with smaller distances between the balls.

Focus on the first ball with both eyes. The goal is to see two incoming and two exiting strings exactly in the centre of the hole of the ball being focussed on. (See illustration). Now gradually focus on the balls behind it. The angle of the string to the tip of the nose and the distances of the balls can be varied as desired to increase the level of difficulty.

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
Pinhole Glasses relax your eyes by reducing the stress factor of light, but the incident light is bundled in a more focussed manner (selectively). This can be used in various ways. In the case of visual stress, which often occurs, for example, in the case of old head injuries or intensive computer work, the glasses can help to improve performance. In addition, Pinhole Glasses provide a clearer view. They are thus often used in combination with other visual exercises. The glasses can be worn 5-10 minutes before the visual relaxation training but also from time to time in everyday life or actively during sports.

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<https://www.youtube.com/@arztz>

Vestibular system

Our balance organs are located in the inner ear. They measure the acceleration of movements of the head and the body and allow us to orient ourselves in space. This information is sent to the brain and forms the basis for our posture and alignment. A functioning balance system is thus the absolute basis for precise execution of movements. The balance system supports the visual system by stabilising the eyes and the visual system supports balance by providing images and thus external fixation points.

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
Basic Exercises Vestibulo-ocular Reflex (VOR)

Train focussing on an object while moving your head asynchronously. To do so, fix the VOR Chart on the wall at eye level and position yourself in the middle with approx. 1 m distance in front of it. Focus your eyes on the little „z“ in the middle of the chart.

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
While your eyes continue to focus on the letter, move your head along the horizontal line to the right and back to the centre (5-10 times each). Repeat the exercise along all the lines shown. Start slowly and increase the tempo, if appropriate. If one axis is particularly difficult for you, you should train it more often.

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© 2023 Ludwig Artzt GmbH